



1505 Bottle, Reagent, Clear, Narrow Mouth, Flathead Stopper

Code No.	Specifications		Packing Quantity	
	Capacity	O.D. x Height x Neck O.D.	Shelf pack	Case
1505S-BT30	30ml	37x70x14mm		10
1505S-BT60	60ml	45x83x14mm		10
1505S-BT120	120ml	54x110x19mm		10
1505S-BT250	250ml	66x130x19mm		10
1505S-BT500	500ml	82x164x24mm		10
1505S-BT1000	1000ml	100x200x29mm		10

●Material is borosilicate glass. Flathead stopper is interchangeable.



1500 Bottle, Reagent, Clear, Narrow Mouth, Flathead Stopper

Code No.	Specifications		Packing Quantity	
	Capacity	O.D. x Height x Neck O.D.	Shelf pack	Case
1500S-BT2000	2000ml	125x225x33mm		6
1500S-BT3000	3000ml	145x290x35mm		6
1500S-BT5000	5000ml	165x350x40mm		1
1500S-BT10000	10000ml	210x440x45mm		1

●Material is borosilicate glass. Flathead stopper is not interchangeable.

51505 Bottle, Reagent, Amber, Narrow Mouth, Flathead Stopper

Code No.	Specifications		Packing Quantity	
	Capacity	O.D. x Height x Neck O.D.	Shelf pack	Case
51505S-BT30	30ml	37x70x14mm		10
51505S-BT60	60ml	45x83x14mm		10
51505S-BT120	120ml	54x110x19mm		10
51505S-BT250	250ml	66x130x19mm		10
51505S-BT500	500ml	82x164x24mm		10
51505S-BT1000	1000ml	100x200x29mm		10

●Material is borosilicate glass. Flathead stopper is interchangeable.



51500 Bottle, Reagent, Amber, Narrow Mouth, Flathead Stopper

Code No.	Specifications		Packing Quantity	
	Capacity	O.D. x Height x Neck O.D.	Shelf pack	Case
51500S-BT2000	2000ml	125x225x33mm		6
51500S-BT3000	3000ml	145x290x35mm		6
51500S-BT5000	5000ml	165x350x40mm		1
51500S-BT10000	10000ml	210x440x45mm		1

●Material is borosilicate glass. Flathead stopper is not interchangeable.



1585 Bottle, Reagent, Clear, Wide Mouth, Flathead Stopper

Code No.	Specifications		Packing Quantity	
	Capacity	O.D. x Height x Neck O.D.	Shelf pack	Case
1585S-BT30	30ml	37x63x24mm		10
1585S-BT60	60ml	45x78x25mm		10
1585S-BT120	120ml	53x100x34mm		10
1585S-BT250	250ml	66x126x45mm		10
1585S-BT500	500ml	82x164x50mm		10
1585S-BT1000	1000ml	100x200x80mm		10

● Material is borosilicate glass. Flathead stopper is interchangeable.

1580 Bottle, Reagent, Clear, Wide Mouth, Flathead Stopper

Code No.	Specifications		Packing Quantity	
	Capacity	O.D. x Height x Neck O.D.	Shelf pack	Case
1580S-BT2000	2000ml	125x250x80mm		6
1580S-BT3000	3000ml	142x285x85mm		6
1580S-BT5000	5000ml	165x340x100mm		1
1580S-BT10000	10000ml	210x440x115mm		1

● Material is borosilicate glass. Flathead stopper is not interchangeable.



51585 Bottle, Reagent, Amber, Wide Mouth, Flathead Stopper

Code No.	Specifications		Packing Quantity	
	Capacity	O.D. x Height x Neck O.D.	Shelf pack	Case
51585S-BT30	30ml	37x63x24mm		10
51585S-BT60	60ml	45x78x25mm		10
51585S-BT120	120ml	53x100x34mm		10
51585S-BT250	250ml	66x126x45mm		10
51585S-BT500	500ml	82x164x50mm		10
51585S-BT1000	1000ml	100x200x80mm		10

● Material is borosilicate glass. Flathead stopper is interchangeable.

51580S Bottle, Reagent, Amber, Wide Mouth, Flathead Stopper

Code No.	Specifications		Packing Quantity	
	Capacity	O.D. x Height x Neck O.D.	Shelf pack	Case
51580S-BT2000	2000ml	125x225x33mm		6
51580S-BT3000	3000ml	145x290x35mm		6
51580S-BT5000	5000ml	165x350x40mm		1
51580S-BT10000	10000ml	210x440x45mm		1

● Material is borosilicate glass. Flathead stopper is not interchangeable.



1700 Bottle, Media, Clear, Narrow Mouth, Screw Cap

Code No.	Capacity	Specifications		Specifications	
		O.D. x Height	Cap size	Shelf pack	Case
1700MEDBT50	50ml	45x100mm	9998CAP430-33		50
1700MEDBT100	100ml	55x129mm	9998CAP430-33		50
1700MEDBT200	200ml	66x148mm	9998CAP430-33		50
1700MEDBT300	300ml	74x164mm	9998CAP430-33		30
1700MEDBT500	500ml	88x192mm	9998CAP430-38		30
1700MEDBT500-33	500ml	88x192mm	9998CAP430-33		30
1700MEDBT1000	1000ml	102x243mm	9998CAP430-38		12

- Screw cap is made of special phenol resin and can withstand repeated sterilization in an autoclave.
- Screw cap is equipped with Teflon® liner which excels in chemical resistance.



51700 Bottle, Media, Amber, Narrow Mouth, Screw Cap

Code No.	Capacity	Specifications		Packing Quantity	
		O.D. x Height	Cap size	Shelf pack	Case
51700MEDBT50	50ml	45x100mm	9998CAP430-33		50
51700MEDBT100	100ml	55x129mm	9998CAP430-33		50
51700MEDBT200	200ml	66x148mm	9998CAP430-33		50
51700MEDBT300	300ml	74x164mm	9998CAP430-33		30
51700MEDBT500	500ml	88x192mm	9998CAP430-38		30
51700MEDBT1000	1000ml	102x243mm	9998CAP430-38		12

- Screw cap is made of special phenol resin and can withstand repeated sterilization in an autoclave.
- Screw cap is equipped with Teflon® liner which excels in chemical resistance.



1705 Bottle, Media, Clear, Narrow Mouth, Screw Cap & Drip-free Pouring Ring

Code No.	Capacity	Specifications		Packing Quantity	
		O.D. x Height	Cap size	Shelf pack	Case
1705MEDBT100	100ml	55x129mm	9998CAP430-33		50
1705MEDBT200	200ml	66x148mm	9998CAP430-33		50
1705MEDBT300	300ml	74x164mm	9998CAP430-33		30
1705MEDBT500	500ml	88x192mm	9998CAP430-38		30
1705MEDBT1000	1000ml	102x243mm	9998CAP430-38		12

- Screw cap is made of special phenol resin and can withstand repeated sterilization in an autoclave.
- Screw cap is equipped with Teflon® liner which excels in chemical resistance.



WARNING

When sterilizing bottles by steam or dry heat, the screw cap must either be removed or just set on the top of the bottle **WITHOUT ENGAGING THE THREAD**. Failure to follow this procedure can result in pressure differences and fracture the bottle.



Bottles



1516 Bottle, Media, Clear, Wide Mouth, Blue Cap

Code No.	Capacity	Specifications		Packing Quantity	
		O.D.x Height x I.D.	STD	Shelf pack	Case
1516BT25	25ml	36x70x14mm	GL-25		10
1516BT50	50ml	46x87x17mm	GL-32		10
1516BT100	100ml	56x100x30mm	GL-45		10
1516BT250	250ml	70x138x30mm	GL-45		10
1516BT500	500ml	86x176x30mm	GL-45		10
1516BT1000	1000ml	101x225x30mm	GL-45		10
1516BT2000	2000ml	136x260x30mm	GL-45		10
1516BT5000	5000ml	181x335x30mm	GL-45		1
1516BT10000	10000ml	234x410x30mm	GL-45		1

- Bottle is made of borosilicate glass and has high thermal resistance and chemical durability.
- Bottle cap and pouring ring can be sterilized up to 140°C.



S1516 Bottle, Media, Clear, Wide Mouth, Red Cap

Code No.	Capacity	Specifications		Packing Quantity	
		O.D.x Height x I.D.	STD	Shelf pack	Case
S1516BT-RED25	25ml	36x70x14mm	GL-25		10
S1516BT-RED50	50ml	46x87x17mm	GL-32		10
S1516BT-RED100	100ml	56x100x30mm	GL-45		10
S1516BT-RED250	250ml	70x138x30mm	GL-45		10
S1516BT-RED500	500ml	86x176x30mm	GL-45		10
S1516BT-RED1000	1000ml	101x225x30mm	GL-45		10
S1516BT-RED2000	2000ml	136x260x30mm	GL-45		10
S1516BT-RED5000	5000ml	181x335x30mm	GL-45		1
S1516BT-RED10000	10000ml	234x410x30mm	GL-45		1

- Bottle is made of borosilicate glass and has high thermal resistance and chemical durability.
- Bottle cap and pouring ring are available for dry heat sterilization (180°C).



1519 Bottle, Media, Amber, Wide Mouth, Blue Cap

Code No.	Capacity	Specifications		Packing Quantity	
		O.D.x Height x I.D.	STD	Shelf pack	Case
1519BT25	25ml	36x70x14mm	GL-25		10
1519BT50	50ml	46x87x17mm	GL-32		10
1519BT100	100ml	56x100x30mm	GL-45		10
1519BT250	250ml	70x138x30mm	GL-45		10
1519BT500	500ml	86x176x30mm	GL-45		10
1519BT1000	1000ml	101x225x30mm	GL-45		10
1519BT2000	2000ml	136x260x30mm	GL-45		10
1519BT5000	5000ml	181x335x30mm	GL-45		1
1519BT10000	10000ml	234x410x30mm	GL-45		1

- Bottle is made of borosilicate glass and has high thermal resistance and chemical durability.
- Bottle cap and pouring ring can be sterilized up to 140°C.
- The bottles are ideal for storage of light sensitive materials.

CAUTION ; Do not autoclave bottles with cap tightly screwed.



S1519 Bottle, Media, Amber, Wide Mouth, Red Cap

Code No.	Capacity	Specifications		Packing Quantity	
		O.D.x Height x I.D.	STD	Shelf pack	Case
S1519BT-RED25	25ml	36x70x14mm	GL-25		10
S1519BT-RED50	50ml	46x87x17mm	GL-32		10
S1519BT-RED100	100ml	56x100x30mm	GL-45		10
S1519BT-RED250	250ml	70x138x30mm	GL-45		10
S1519BT-RED500	500ml	86x176x30mm	GL-45		10
S1519BT-RED1000	1000ml	101x225x30mm	GL-45		10
S1519BT-RED2000	2000ml	136x260x30mm	GL-45		10
S1519BT-RED5000	5000ml	181x335x30mm	GL-45		1
S1519BT-RED10000	10000ml	234x410x30mm	GL-45		1

- Bottle is made of borosilicate glass and has high thermal resistance and chemical durability.
- Bottle cap and pouring ring are available for dry heat sterilization (180°C).
- The bottles are ideal for storage of light sensitive materials.



1518 Bottle, Media, Wide Mouth, Safety Coated

Code No.	Capacity	Specifications		Packing Quantity	
		O.D.x Height x I.D.	STD	Shelf pack	Case
1518BT100	100ml	56x100x30mm	GL-45		10
1518BT250	250ml	70x138x30mm	GL-45		10
1518BT500	500ml	86x176x30mm	GL-45		10
1518BT1000	1000ml	101x225x30mm	GL-45		10
1518BT2000	2000ml	136x260x30mm	GL-45		10
1518BT5000	5000ml	181x335x30mm	GL-45		1

- Bottle is made of borosilicate glass and has high thermal resistance and chemical durability.
- Plastic coating provides protection against mechanical breakage. The Co-polymer based coating is bonded firmly to the glass and remains transparent even after repeated steam sterilization at upto 121°C.



1570 Bottle, Media, Wide Mouth, Square Shape

Code No.	Capacity	Specifications		Packing Quantity	
		O.D.x Height	STD	Shelf pack	Case
1570BT100	100ml	30x109mm	GL-32		10
1570BT250	250ml	64x143mm	GL-45		10
1570BT500	500ml	78x181mm	GL-45		10
1570BT1000	1000ml	94x222mm	GL-45		10

- Square shape allows closer stacking in water baths or on shelves
- Up to 20% saving of space compared to round bottles
- Bottle complete with cap and pouring ring
- Withstands steam sterilization or dry heat up to 140°C
- Wide neck for easy filling



WARNING

When sterilizing bottles by steam or dry heat, the screw cap must either be removed or just set on the top of the bottle **WITHOUT ENGAGING THE THREAD**. Failure to follow this procedure can result in pressure differences and fracture the bottle.

CAUTION ; Do not autoclave bottles with cap tightly screwed.



Bottles



4506 Screw Cap for 1516 Media Bottle

Code No.	Specifications		Packing Quantity
	Thread size	Material	
4506CAP25D	GL-25	Polypropylene	10
4506CAP32D	GL-32	Polypropylene	10
4506CAP45D	GL-45	Polypropylene	10

● Color; Light blue. Autoclavable up to 140°C



4526 Screw Cap for S1516 Media Bottle

Code No.	Specifications		Packing Quantity
	Thread size	Material	
4526CAP25D	GL-25	Thermoplastic polyester	10
4526CAP32D	GL-32	Thermoplastic polyester	10
4526CAP45D	GL-45	Thermoplastic polyester	10

● Color ; Red. Can withstand up to 180°C hot air sterilization. Complete with PTFE protective seal.

Note:

Color of cap may discolor when used at maximum temperature of 180°C. However, such discoloration does not affect performance.



4506 Colored Cap for Media bottle, Wide Mouth

Code No.	Specifications			Packing Quantity
	Thread size	Material	Color	
4506CAP45E	GL-45	Polypropylene	Green	10
4506CAP45F	GL-45	Polypropylene	Yellow	10
4506CAP45G	GL-45	Polypropylene	Gray	10

● Can withstand up to 140°C



4516 Spare Pouring Ring

Code No.	Specifications		Packing Quantity
	Thread size	Material	
4516RING32D	GL-32	Polypropylene	10
4516RING45D	GL-45	Polypropylene	10

● Color; Clear. Can withstand up to 140°C



4536 Spare Pouring Ring

Code No.	Specifications		Packing Quantity
	Thread size	Material	
4536RING32D	GL-32	PTFE	10
4536RING45D	GL-45	PTFE	10

● Color; Red. Can withstand up to 180°C



1650 Bottle, Incubation

Code No.	Specification		Packing Quantity	
	Capacity	Stopper (Standard Taper)	Shelf pack	Case
1650FBT100	100ml	19		25
1650FBT200	200ml	19		10
1650FBT300	300ml	19		10



1651 Bottle, Incubation, Glass Cup

Code No.	Specification		Packing Quantity	
	Capacity	Stopper (Standard Taper)	Shelf pack	Case
1651FBT100	100ml	19		25
1651FBT200	200ml	19		10
1651FBT300	300ml	19		10



1654 Bottle, Incubation, Rubber Cup

Code No.	Specification		Packing Quantity	
	Capacity	Stopper (Standard Taper)	Shelf pack	Case
1654FBT100	100ml	19		10
1654FBT200	200ml	24		10
1654FBT300	300ml	24		10



Bottle, Weighing

Code No.	Specification	Packing Quantity	
	O.D. x Height	Shelf pack	Case
WEIGHBT18-30	18x30mm		10
WEIGHBT25-45	25x45mm		10
WEIGHBT30-30	30x30mm		10
WEIGHBT30-45	30x45mm		10
WEIGHBT30-60	30x60mm		10
WEIGHBT40-40	40x40mm		10
WEIGHBT40-60	40x60mm		10
WEIGHBT45-45	45x45mm		10
WEIGHBT50-30	50x30mm		10
WEIGHBT50-50	50x50mm		10
WEIGHBT50-75	50x75mm		10
WEIGHBT60-40	60x40mm		10
WEIGHBT60-75	60x75mm		10



3161 Dish, Culture, Petri (Sheet Glass made)

Code No	Size	Specifications	Packing Quantity	
		O.D. x Height (Approx.)	Shelf Pack	Case
3161DISH90-20	90-20	97x17mm		50

