



## 3081 Desiccator, Knob Top, Porcelain Plate

Code No	Specifications (Approx.)			Packing Quantity	
	Body I.D.	Overall Height	Disc Diameter	Shelf pack	Case
3081 DSC120	120mm	210mm	105mm		1
3081 DSC150	150mm	230mm	135mm		1
3081 DSC180	180mm	280mm	165mm		1
3081 DSC210	210mm	330mm	195mm		1
3081 DSC240	240mm	340mm	225mm		1
3081 DSC300	300mm	430mm	285mm		1

- Not available for vacuum applications.
- Do not heat or subject to pressure.



## 3121 Desiccator, Opening for Rubber Stopper, Porcelain Plate

Code No	Specifications (Approx.)			Packing Quantity	
	Body I.D.	Overall Height	Disc Diameter	Shelf pack	Case
3121 DSC120	120mm	210mm	105mm		1
3121 DSC150	150mm	230mm	135mm		1
3121 DSC180	180mm	280mm	165mm		1
3121 DSC210	210mm	330mm	195mm		1
3121 DSC240	240mm	340mm	225mm		1
3121 DSC300	300mm	430mm	285mm		1

- For sustained vacuum storage, all ground surfaces should be cleaned thoroughly and pasted a silicon based-vacuum grease.
- Do not heat or subject to pressure.



## Stopcock for 3121 Desiccator, Rubber Stopper

Code No	Specifications (Approx.)			Packing Quantity	
	Rubber stopper	Nozzle O.D.	Stopcock	Shelf pack	Case
3121 SDOCK7	NO. 10 (for 120,150)	11 mm	NO. 16		1
3121 SDOCK9	NO. 11 (for 180)	11 mm	NO. 16		1
3121 SDOCK10	NO. 13 (for 210,240)	11 mm	NO. 16		1
3121 SDOCK13	NO. 14 (for 300)	11 mm	NO. 16		1

- O.D. of glass tubing ; 9mm



## Desiccator Plate

Code No.	Specifications		Packing Quantity	
	Material	Plate O.D.	Shelf pack	Case
DSC-P120G	Porcelain	105mm		1
DSC-P150G	Porcelain	135mm		1
DSC-P180G	Porcelain	165mm		1
DSC-P210G	Porcelain	195mm		1
DSC-P240G	Porcelain	225mm		1
DSC-P300G	Porcelain	285mm		1
DSC-P180S	Stainless Steel	165mm		1
DSC-P210S	Stainless Steel	195mm		1
DSC-P240S	Stainless Steel	225mm		1

