



## Dispenser, Clear, Amber Bottle

Code No.	Specifications			Packing Quantity		
	Capacity	Graduation	Tolerance	Bottle	Shelf pack	Case
S-BUNCHU-SET1	1ml	0.01ml	±1%	1000ml		1
S-BUNCHU-SET2	2ml	0.02ml	±1%	1000ml		1
S-BUNCHU-SET5	5ml	0.1ml	±1%	1000ml		1
S-BUNCHU-SET10	10ml	0.1ml	±1%	1000ml		1
S-BUNCHU-SET20	20ml	0.2ml	±1%	1000ml		1
S-BUNCHU-SET25	25ml	0.2ml	±1%	1000ml		1
S-BUNCHU-SET50	50ml	1.0ml	±1%	3800ml		1
S-BUNCHU-SET100	100ml	1.0ml	±1%	3800ml		1

- Visible blue color graduations.
- Heavy weight bottle for suitable operation.



## Dispenser, Amber, with Bottle

Code No.	Specifications			Packing Quantity		
	Capacity	Graduation	Tolerance	Bottle	Shelf pack	Case
S-B-BUNCHU-SET1	1ml	0.01ml	±1%	1000ml		1
S-B-BUNCHU-SET2	2ml	0.02ml	±1%	1000ml		1
S-B-BUNCHU-SET5	5ml	0.1ml	±1%	1000ml		1
S-B-BUNCHU-SET10	10ml	0.1ml	±1%	1000ml		1
S-B-BUNCHU-SET20	20ml	0.2ml	±1%	1000ml		1
S-B-BUNCHU-SET25	25ml	0.2ml	±1%	1000ml		1
S-B-BUNCHU-SET50	50ml	1.0ml	±1%	3800ml		1
S-B-BUNCHU-SET100	100ml	1.0ml	±1%	3800ml		1



## Dispenser, Clear, Silicon Rubber, without Bottle

Code No.	Specifications				Packing Quantity	
	Capacity	Graduation	Tolerance	Silicon Rubber	Shelf pack	Case
PS-BUNCHU1S	1ml	0.01ml	±1%	No.4		1
PS-BUNCHU2S	2ml	0.02ml	±1%	No.4		1
PS-BUNCHU5S	5ml	0.1ml	±1%	No.4		1
PS-BUNCHU10S	10ml	0.1ml	±1%	No.4		1
PS-BUNCHU20S	20ml	0.2ml	±1%	No.8		1
PS-BUNCHU25S	25ml	0.2ml	±1%	No.8		1
PS-BUNCHU50S	50ml	1.0ml	±1%	No.15		1
PS-BUNCHU100S	100ml	1.0ml	±1%	No.15		1

- Available for use with erlenmyer flask, gallon bottle, etc.

## Dispenser, Amber, Silicon Rubber, without Bottle

Code No.	Specifications				Packing Quantity	
	Capacity	Graduation	Tolerance	Silicon Rubber	Shelf pack	Case
PS-B-BUNCHU1S	1ml	0.01ml	±1%	No.4		1
PS-B-BUNCHU2S	2ml	0.02ml	±1%	No.4		1
PS-B-BUNCHU5S	5ml	0.1ml	±1%	No.4		1
PS-B-BUNCHU10S	10ml	0.1ml	±1%	No.4		1
PS-B-BUNCHU20S	20ml	0.2ml	±1%	No.8		1
PS-B-BUNCHU25S	25ml	0.2ml	±1%	No.8		1
PS-B-BUNCHU50S	50ml	1.0ml	±1%	No.15		1
PS-B-BUNCHU100S	100ml	1.0ml	±1%	No.15		1

- Available for use with erlenmyer flask, gallon bottle, etc.





# Dispenser



## Bottle, Amber, for Dispenser

Code No.	Description	Packing Quantity	
		Shelf pack	Case
S-BUN-BT1000	1000ml		1
S-BUN-BT2000	2000ml		1
S-BUN-BT3800	3800ml		1



## 1820 Dispenser, Amber Bottle

Code No.	Specifications				Packing Quantity	
	Capacity	Graduation	Tolerance	Bottle	Shelf pack	Case
1820BUNCHU0.5	0.5ml	0.01ml	±1%	500ml		1
1820BUNCHU1	1ml	0.02ml	±1%	500ml		1
1820BUNCHU2	2ml	0.05ml	±1%	500ml		1
1820BUNCHU5	5ml	0.1ml	±1%	500ml		1
1820BUNCHU10	10ml	0.2ml	±1%	1000ml		1

- All portions in direct contact with liquid are made of glass and PTFE.
- Injection cylinder is included in the bottle to prevent breakage.
- Liquid replenishment into bottle can be easily made with funnel attachment without opening the cap.