


Flasks with joint, Stopper




4323A Flask, Boiling, Round Bottom, Short Neck, Two Necks

| Code No | Specifications | | | Packing Quantity | |
|--------------|----------------|--------------|------------|------------------|------|
| | Capacity | Center joint | Side joint | Shelf pack | Case |
| 4323A-100-24 | 100ml | 24/29 | 14/23 | | 1 |
| 4323A-100-29 | 100ml | 29/32 | 14/23 | | 1 |
| 4323A-250-24 | 250ml | 24/29 | 14/23 | | 1 |
| 4323A-250-29 | 250ml | 29/32 | 14/23 | | 1 |
| 4323A-500-24 | 500ml | 24/29 | 14/23 | | 1 |
| 4323A-500-29 | 500ml | 29/32 | 14/23 | | 1 |
| 4323A-1L-24 | 1000ml | 24/29 | 14/23 | | 1 |
| 4323A-1L-29 | 1000ml | 29/32 | 14/23 | | 1 |
| 4323A-2L-29 | 2000ml | 29/32 | 14/23 | | 1 |



4323B Flask, Boiling, Round Bottom, Short Neck, Three Necks

| Code No | Specifications | | | Packing Quantity | |
|--------------|----------------|--------------|------------|------------------|------|
| | Capacity | Center joint | Side joint | Shelf pack | Case |
| 4323B-100-29 | 100ml | 29/32 | 14/23 | | 1 |
| 4323B-250-24 | 250ml | 24/29 | 14/23 | | 1 |
| 4323B-500-24 | 500ml | 24/29 | 14/23 | | 1 |
| 4323B-500-29 | 500ml | 29/32 | 14/23 | | 1 |
| 4323B-1L-24 | 1000ml | 24/29 | 14/23 | | 1 |
| 4323B-1L-29 | 1000ml | 29/32 | 14/23 | | 1 |



7570F Stopper, Flat top

| Code No | Specifications | Packing Quantity | |
|--------------------|----------------|------------------|------|
| | Joint Size | Shelf pack | Case |
| 7570F-STOPER 15/25 | 15/25 | | 100 |
| 7570F-STOPER 15/35 | 15/35 | | 100 |
| 7570F-STOPER 19/28 | 19/28 | | 100 |
| 7570F-STOPER 19/38 | 19/38 | | 100 |
| 7570F-STOPER 24/30 | 24/30 | | 100 |
| 7570F-STOPER 24/40 | 24/40 | | 100 |
| 7570F-STOPER 29/32 | 29/32 | | 50 |
| 7570F-STOPER 29/42 | 29/42 | | 50 |
| 7570F-STOPER 9 | 9 | | 100 |
| 7570F-STOPER 13 | 13 | | 100 |
| 7570F-STOPER 16 | 16 | | 100 |
| 7570F-STOPER 19 | 19 | | 100 |
| 7570F-STOPER 22 | 22 | | 100 |
| 7570F-STOPER 27 | 27 | | 50 |
| 7570F-STOPER 38 | 38 | | 50 |

