


Joins (J.I.S. Standard)



6560JT Ground Joint, Full Length, Inner Part, Interchangeable

Code No.	Specifications				Packing Quantity	
	Joint	Joint O.D.	Joint Length	Tube O.D.	Shelf pack	Case
6560JT-I 10/30	10/30	10.0mm	30mm	8mm		10
6560JT-I 12/30	12/30	12.0mm	30mm	10mm		10
6560JT-I 15/25	15/25	15.5mm	25mm	13mm		10
6560JT-I 15/35	15/35	15.5mm	35mm	13mm		10
6560JT-I 19/28-16	19/28	18.8mm	28mm	16mm		10
6560JT-I 19/38-16	19/38	18.8mm	38mm	16mm		10
6560JT-I 24/30-20	24/30	24.0mm	30mm	20mm		10
6560JT-I 24/40-20	24/40	24.0mm	40mm	20mm		10
6560JT-I 29/32-25	29/32	29.2mm	32mm	25mm		10
6560JT-I 29/42-25	29/42	29.2mm	42mm	25mm		10
6560JT-I 34/45-30	34/45	34.5mm	45mm	30mm		10
6560JT-I 40/50	40/50	40.0mm	50mm	38mm		10
6560JT-I 45/50	45/50	45.0mm	50mm	40mm		10

●J.I.S. ; Abbreviation of Japan Industrial Standard ,nearly same as ASTM specifications.



6580JT Ground Joint, Full Length, Outer Part, Interchangeable

Code No.	Specifications				Packing Quantity	
	Joint	Joint O.D.	Joint Length	Tube O.D.	Shelf pack	Case
6580JT-O 10/30	10/30	10.0mm	30mm	14mm		10
6580JT-O 12/30	12/30	12.0mm	30mm	16mm		10
6580JT-O 15/25	15/25	15.5mm	25mm	18mm		10
6580JT-O 15/35	15/35	15.5mm	35mm	18mm		10
6580JT-O 19/28	19/28	18.8mm	28mm	22mm		10
6580JT-O 19/38	19/38	18.8mm	38mm	22mm		10
6580JT-O 24/30	24/30	24.0mm	30mm	28mm		10
6580JT-O 24/40	24/40	24.0mm	40mm	28mm		10
6580JT-O 29/32	29/32	29.2mm	32mm	32mm		10
6580JT-O 29/42	29/42	29.2mm	42mm	32mm		10
6580JT-O 34/45	34/45	34.5mm	45mm	38mm		10
6580JT-O 40/50	40/50	40.0mm	50mm	45mm		10
6580JT-O 45/50	45/50	45.0mm	50mm	48mm		10

●J.I.S. ; Abbreviation of Japan Industrial Standard ,nearly same as ASTM specifications.

6762GJ Ball Joint

Code No.	Specifications				Packing Quantity	
	Joint	Joint O.D.	Hole Diameter	Tube O.D.	Shelf pack	Case
6762GJ-I 13/1	S13/1	12.7mm	1mm	8mm		10
6762GJ-I 13/2	S13/2	12.7mm	2mm	8mm		10
6762GJ-I 13/3	S13/3	12.7mm	3mm	9mm		10
6762GJ-I 13/5	S13/5	12.7mm	5mm	9mm		10
6762GJ-I 19/7	S19/7	19.05mm	7mm	10mm		10
6762GJ-I 19/9	S19/9	19.05mm	9mm	12mm		10
6762GJ-I 22/5	S22/5	22.225mm	5mm	10mm		10
6762GJ-I 22//9	S22/9	22.225mm	9mm	13mm		10
6762GJ-I 29/12	S29/12	28.575mm	12mm	16mm		10
6762GJ-I 29/15	S29/15	28.575mm	15mm	20mm		10
6762GJ-I 35/20	S35/20	34.925mm	20mm	25mm		10

●J.I.S. ; Abbreviation of Japan Industrial Standard ,nearly same as ASTM specifications.

6764GJ Socket Joint

Code No.	Specifications				Packing Quantity	
	Joint	Joint O.D.	Hole Diameter	Tube O.D.	Shelf pack	Case
6764GJ-O 13/1	S13/1	12.7mm	1mm	8mm		10
6764GJ-O 13/2	S13/2	12.7mm	2mm	8mm		10
6764GJ-O 13/3	S13/3	12.7mm	3mm	9mm		10
6764GJ-O 13/5	S13/5	12.7mm	5mm	9mm		10
6764GJ-O 19/7	S19/7	19.05mm	7mm	10mm		10
6764GJ-O 19/9	S19/9	19.05mm	9mm	12mm		10
6764GJ-O 22/5	S22/5	22.225mm	5mm	10mm		10
6764GJ-O 22//9	S22/9	22.225mm	9mm	13mm		10
6764GJ-O 29/12	S29/12	28.575mm	12mm	16mm		10
6764GJ-O 29/15	S29/15	28.575mm	15mm	20mm		10
6764GJ-O 35/20	S35/20	34.925mm	20mm	25mm		10

●J.I.S. ; Abbreviation of Japan Industrial Standard ,nearly same as ASTM specifications.

