




# Joints, Adapter-I.S.O. Standard (Made to order)



## 6560 Ground Joint, Full Length, Inner Part, Interchangeable


Code No	Specifications			Packing Quantity	
	Joint Size	Ground Portion O.D	Ground Portion Length	Shelf pack	Case
6560-7	7/16	7.5mm	16mm		10
6560-10	10/19	10mm	19mm		10
6560-12	12/21	12.5mm	21mm		10
6560-14	14/23	14.5mm	23mm		10
6560-19	19/26	18.8mm	26mm		10
6560-24	24/29	24mm	29mm		10
6560-29	29/32	29.2mm	32mm		10
6560-34	34/45	34.5mm	35mm		10
6560-45	45/40	45mm	40mm		10

made to order

## 6580 Ground Joint, Full Length, Outer Part,

Code No	Specifications			Packing Quantity	
	Joint Size	Ground Portion O.D	Ground Portion Length	Shelf pack	Case
6580-7	7/16	7.5mm	16mm		10
6580-10	10/19	10mm	19mm		10
6580-12	12/21	12.5mm	21mm		10
6580-14	14/23	14.5mm	23mm		10
6580-19	19/26	18.8mm	26mm		10
6580-24	24/29	24mm	29mm		10
6580-29	29/32	29.2mm	32mm		10
6580-34	34/45	34.5mm	35mm		10
6580-45	45/40	45mm	40mm		10

made to order



## 8800 Tube, Adapter, Reducing

Code No	Specifications			Packing Quantity	
	Size	Outer Joint	Inner joint	Shelf pack	Case
8800-14-19	14-19-23	14/23	19/26		1
8800-14-24	14-24-23	14/23	24/29		1
8800-14-29	14-29-23	14/23	29/32		1
8800-19-14	19-14-26	19/26	14/23		1
8800-19-24	19-24-26	19/26	24/29		1
8800-19-29	19-29-26	19/26	29/32		1
8800-24-14	24-14-29	24/29	14/23		1
8800-24-19	24-19-29	24/29	29/26		1
8800-24-29	24-29-29	24/29	29/32		1
8800-29-14	29-14-32	29/32	14/23		1
8800-29-19	29-19-32	29/32	19/26		1
8800-29-24	29-24-32	29/32	24/29		1

made to order

## 8820 Tube, Adapter, Enlarging

Code No	Specifications			Packing Quantity	
	Size	Outer joint	Inner joint	Shelf pack	Case
8820-14-19	14-19-23	19/26	14/23		1
8820-14-24	14-24-23	24/29	14/23		1
8820-14-29	14-29-23	29/32	14/23		1
8820-19-14	19-14-26	14/23	19/26		1
8820-19-24	19-24-26	24/29	19/26		1
8820-19-29	19-29-26	29/32	19/26		1
8820-24-14	24-14-29	14/23	24/29		1
8820-24-19	24-19-29	29/26	24/29		1
8820-24-29	24-29-29	29/32	24/29		1
8820-29-14	29-14-32	14/23	29/32		1
8820-29-19	29-19-32	19/26	29/32		1
8820-29-24	29-24-32	24/29	29/32		1

made to order