



Plant Tissue Culture Ware



Flask, for Meristem Culture

Code No	Specifications		Packing Quantity	
	Capacity	Widthx Heightx Neck I.D.	Shelf pack	Case
MER-FK300	300ml	90x127x40mm		60
MER-FK500	500ml	108x159x40mm		48

- These flasks have wide mouths and are ideal for Meristem culture.
- Wide mouth facilitates transfer of plants. All flasks have permanent white enamel graduations and marking spots for easy identification.



Flask, For Plant Tissue Culture

Code No	Specifications		Packing Quantity	
	Capacity	Widthx Heightx Neck I.D.	Shelf pack	Case
CUL-FK100	100ml	59x90x47mm		90
CUL-FK200	200ml	74x110x54mm		60
CUL-FK300	300ml	84x135x58mm		56

- These flasks have wide mouths and are ideal for plant tissue culture.
- All flasks have permanent white enamel graduations and marking spots for easy identification.



4442 Flask, for Plant Tissue Culture, Long Neck

Code No	Specifications		Packing Quantity	
	Capacity	Widthx Heightx Neck I.D.	Shelf pack	Case
4442FK100	100ml	65x125x25mm		12
4442FK200	200ml	81x146x38mm		12
4442FK300	300ml	91x165x38mm		12
4442FK500	500ml	108x188x38mm		12

- This long neck flask is erlenmeyer shape to which a long straight and rimless neck has been attached.
- The longer neck reduces amount of splashed solution reaching the closure. It is ideal for shaker cultures.



Plant Box, Plastic Jar for Plant Tissue Culture

Code No	Specifications			Packing Quantity	
	Capacity	Widthx Depthx Height	Shelf pack	Case	
CUL-JAR300	300ml	60x60x100mm		100	

- This reusable plant box is made of polycarbonate and can withstand repeated sterilization in an autoclave.



Test Tube, with Rim, Flat Bottom

Code No	Specifications		Packing Quantity	
	O.D. x Length	Thickness	Shelf pack	Case
TEST-F25-100	25x100mm	1.2mm		100
TEST-F25-120	25x120mm	1.2mm		100
TEST-F25-150	25x150mm	1.2mm		100
TEST-F30-120	30x120mm	1.2mm		100
TEST-F40-130	40x130mm	2.0mm		100



Test Tube, without Rim, Flat Bottom

Code No	Specifications		Packing Quantity	
	O.D. x Length	Thickness	Shelf pack	Case
9820TST-F25-100	25x100mm	1.2mm		100
9820TST-F25-120	25x120mm	1.2mm		100
9820TST-F25-150	25x150mm	1.2mm		100
9820TST-F30-120	30x120mm	1.2mm		100
9820TST-F40-130	40x130mm	2.0mm		100



Cap for Test Tube without Rim and 4442 Flask

Code No	Description	Packing Quantity	
		Shelf pack	Case
9988CAP25	For 9820TST-F25	100	1,000
9988CAP30	For 9820TST-F30	100	1,000
9988CAP38	For 4442FK		100
9988CAP40	For 9820TST-F40		100

● This cap is made of polypropylene and can withstand sterilization in an autoclave.





Plant Tissue Culture Ware



7601 Stopper, Silicone Rubber

Code No	Specifications				Packing Quantity	
	Top O.D.	Bottom O.D.	Height	Hole Dia.	Shelf pack	Case
7601ST015	43mm	38.5mm	25mm	6mm		100

● This stopper is made of silicone rubber and can withstand repeated sterilization in an autoclave.

7621 Stopper, Rubber

Code No	Specifications				Packing Quantity	
	Top O.D.	Bottom O.D.	Height	Hole Dia.	Shelf pack	Case
7621ST015	43mm	38.5mm	25mm	6mm		100

● This stopper is made of black rubber and can withstand several sterilization in an autoclave.



Test Tube Rack

Code No	Description	Packing Quantity	
		Shelf pack	Case
TEST-RACK25-50	50pcs for 25mm Test Tube	1	10
TEST-RACK30-50	50pcs for 30mm Test Tube	1	10
TEST-RACK40-10	10pcs for 40mm Test Tube	1	10
TEST-RACK40-50	50pcs for 40mm Test Tube	1	10
TEST-RACK25-36	36pcs for 25mm Test Tube	1	10
TEST-RACK30-25	25pcs for 30mm Test Tube	1	10
TEST-RACK40-16	16pcs for 40mm Test Tube	1	10

● This rack is made of stainless steel and can withstand sterilization in an autoclave.



Pipette, Callus Transfer, Large Tip Opening


Code No	Specifications			Packing Quantity	
	Capacity	Length	I.D. at Tip	Shelf pack	Case
K-PIPET-LT2	2ml	245mm	3mm		10
K-PIPET-LT5	5ml	275mm	4mm		10
K-PIPET-LT10	10ml	340mm	4mm		10

● Suitable for transfer of callus suspension media.



Surgical Knife

Code No	Description	Packing Quantity	
		Shelf pack	Case
MES-HOL3	Holder only	1	5
MES-HOL4	Holder only	1	5
MES-BR11	Blade only, for MES-HOL3, sharp		20
MES-BR21	Blade only, for MES-HOL3, round		20



Needle, with holder

Code No	Description	Packing Quantity	
		Shelf pack	Case
NED-HOL-A	Needle with holder, Aluminum	1	10

Forceps

Code No	Specifications		Packing Quantity	
	Description	Length	Shelf pack	Case
PIN1 8-190	Straight	190mm	1	10
PIN1 8-240	Straight	240mm	1	1
PIN1 8-300	Straight	300mm	1	1
PIN2 8-175	Straight	175mm	1	1
PIN2 8-250	Straight	250mm	1	1
PIN2 8-300	Straight	300mm	1	1
PIN1 9-190	Bend	190mm	1	10
PIN1 9-240	Bend	240mm	1	1
PIN1 9-300	Bend	300mm	1	1
PIN1 7-150A	Straight	150mm	1	10
PIN1 7-230A	Straight	230mm	1	10
PIN2 2-125	Straight	125mm	1	10

