



1880 Vial, Sample, Black Phenolic Cap

Code No.	Specification			Packing Quantity	
	Capacity	OD x Height x Neck ID	Cap	Shelf pack	Case
1880SV2	2ml	15x40x7.4mm	9996CAP13	100	600
1880SV3	3ml	15x45x7.4mm	9996CAP13	100	600
1880SV5	5ml	16.5x45x9.0mm	9996CAP15	100	600
1880SV7	7ml	16.5x55x9.0mm	9996CAP15	100	600
1880SV10	10ml	19x55x9.0mm	9996CAP15	100	600
1880SV15	15ml	25x55x13.1mm	9996CAP20	50	200
1880SV20	20ml	27.5x55x17.4mm	9996CAP24	50	200
1880SV30	30ml	27.5x70x17.4mm	9996CAP24	50	200

- Cap is made of special phenol resin and can withstand repeated sterilization in an autoclave.
- Under freezing, sample should be limited to one fourth of the nominal capacity.
- This vial should not be used for the solvent such as chloroform.



1890 Vial, Sample, White Polypropylene Cap


Code No.	Specification			Packing Quantity	
	Capacity	OD x Height x Neck ID	Cap	Shelf pack	Case
1890CV2	2ml	15x40x7.4mm	9995CAP13	100	600
1890CV3	3ml	15x45x7.4mm	9995CAP13	100	600
1890CV5	5ml	16.5x45x9.0mm	9995CAP15	100	600
1890CV7	7ml	16.5x55x9.0mm	9995CAP15	100	600
1890CV10	10ml	19x55x9.0mm	9995CAP15	100	600
1890CV15	15ml	25x55x13.1mm	9995CAP20	50	200
1890CV20	20ml	27.5x55x17.4mm	9995CAP24	50	200
1890CV30	30ml	27.5x70x17.4mm	9995CAP24	50	200
1890CV50	50ml	40x70x20.1mm	9995CAP33	50	200
1890CV100	100ml	40x120x20.1mm	9995CAP33	50	200

- Polypropylene cap with polypropylene liner provides excellent resistance against acid, alkali and other solvents.
- Under freezing, sample should be limited to one fourth of the nominal capacity.



1870 Vial, Snap Cap

Code No.	Specification			Packing Quantity	
	Capacity	OD x Height x Neck ID	Cap	Shelf pack	Case
1870SNV3	3ml	16.5x40x12.8mm	9994CAP2	100	900
1870SNV5	5ml	18x40x12.8mm	9994CAP2	100	900
1870SNV10	10ml	21x45x14.5mm	9994CAP3	100	600
1870SNV15	15ml	24x50x17.0mm	9994CAP4	50	400
1870SNV20	20ml	27.5x55x17.0mm	9994CAP4	50	400
1870SNV30	30ml	30x63x20.7mm	9994CAP5	50	400
1870SNV50	50ml	40x75x25.5mm	9994CAP7	50	200
1870SNV100	100ml	40x120x25.5mm	9994CAP7	50	200



1855 Vial, Scintillation, Screw Cap

Code No.	Specification			Packing Quantity	
	Capacity	OD x Height	Cap	Shelf pack	Case
1855SIV7	7ml	16.5x54mm	9997CAP425-15	100	500
1855SIV20	20ml	27.5x57mm	9997CAP400-22	100	500

- Five inner boxes are packed in one outer carton.
- Inner box without cover provides easy handling and vial is capped to avoid contamination.

