




Volumetric Ware Buret




Buret, Straight Bore, Glass Stopcock, Class A

Code No.	Specifications				Packing Quantity	
	Capacity	Graduation	Tolerance	Stopcock	Shelf pack	Case
BURET10S	10ml	0.05ml	±0.02ml	No.13		1
BURET25S	25ml	0.1ml	±0.03ml	No.13		1
BURET50S	50ml	0.1ml	±0.05ml	No.13		1
BURET100S	100ml	0.2ml	±0.1ml	No.13		1

Buret, Straight Bore, Amber, Glass Stopcock, Class A

Code No.	Specifications				Packing Quantity	
	Capacity	Graduation	Tolerance	Stopcock	Shelf pack	Case
B-BURET10S	10ml	0.05ml	±0.02ml	No.13		1
B-BURET25S	25ml	0.1ml	±0.03ml	No.13		1
B-BURET50S	50ml	0.1ml	±0.05ml	No.13		1
B-BURET100S	100ml	0.2ml	±0.1ml	No.13		1



2103 Buret, Straight Bore, Teflon® Stopcock, Class A

Code No.	Specifications				Packing Quantity	
	Capacity	Graduation	Tolerance	Stopcock	Shelf pack	Case
2103BURET10S	10ml	0.05ml	±0.02ml	No.13.5		1
2103BURET25S	25ml	0.1ml	±0.03ml	No.13.5		1
2103BURET50S	50ml	0.1ml	±0.05ml	No.13.5		1

B-2103 Buret, Straight Bore, Amber, Teflon Stopcock, Class A


Code No.	Specifications				Packing Quantity	
	Capacity	Graduation	Tolerance	Stopcock	Shelf pack	Case
B-2103BURET10S	10ml	0.05ml	±0.02ml	No.13.5		1
B-2103BURET25S	25ml	0.1ml	±0.03ml	No.13.5		1
B-2103BURET50S	50ml	0.1ml	±0.05ml	No.13.5		1



Buret, Automatic, Glass Stopcock and Flat Joint Bottle

Code No.	Specifications				Packing Quantity	
	Capacity	Graduation	Tolerance	Bottle	Shelf pack	Case
F-BURET10S	10ml	0.05ml	±0.02ml	1000ml		1
F-BURET25S	25ml	0.1ml	±0.03ml	2000ml		1
F-BURET50S	50ml	0.1ml	±0.05ml	2000ml		1


- Easy to remove bottle from body.
- Flat joint portion is transparent. Hard to damage and easy to remove stains from bottle.
- All necessary parts for assembly are completely attached.



Buret, Automatic, Amber, Glass Stopcock and Flat Joint Bottle

Code No.	Specifications				Packing Quantity	
	Capacity	Graduation	Tolerance	Bottle	Shelf pack	Case
F-B-BURET10S	10ml	0.05ml	±0.02ml	1000ml		1
F-B-BURET25S	25ml	0.1ml	±0.03ml	2000ml		1
F-B-BURET50S	50ml	0.1ml	±0.05ml	2000ml		1

- Easy to remove bottle from body.
- Flat joint portion is transparent. Hard to damage and easy to remove stains from bottle.
- All necessary parts for assembly are completely attached.



Buret, Automatic, Teflon® Stopcock and Flat Joint Bottle

Code No.	Specifications				Packing Quantity	
	Capacity	Graduation	Tolerance	Bottle	Shelf pack	Case
F-2103BURET10S	10ml	0.05ml	±0.02ml	1000ml		1
F-2103BURET25S	25ml	0.1ml	±0.03ml	2000ml		1
F-2103BURET50S	50ml	0.1ml	±0.05ml	2000ml		1

- Easy to remove bottle from body.
- Flat joint portion is transparent. Hard to damage and easy to remove stains from bottle.
- All necessary parts for assembly are completely attached.

Buret, Automatic, Amber, Teflon® Stopcock and Flat Joint Bottle

Code No.	Specifications				Packing Quantity	
	Capacity	Graduation	Tolerance	Bottle	Shelf pack	Case
F-B-2103BURET10S	10ml	0.05ml	±0.02ml	1000ml		1
F-B-2103BURET25S	25ml	0.1ml	±0.03ml	2000ml		1
F-B-2103BURET50S	50ml	0.1ml	±0.05ml	2000ml		1

- Easy to remove bottle from body.
- Flat joint portion is transparent. Hard to damage and easy to remove stains from bottle.
- All necessary parts for assembly are completely attached.




Bottle for Buret, Automatic, Flat Joint

Code No.	Specification	Packing Quantity	
		Shelf pack	Case
F-BTL1000	Bottle, Clear 1000ml		1
F-BTL2000	Bottle, Clear 2000ml		1
F-B-BTL1000	Bottle, Amber, 1000ml		1
F-B-BTL2000	Bottle, Amber, 2000ml		1



Volumetric Ware Buret



Buret, Micro, without Stand


Code No	Specifications			Packing Quantity	
	Capacity	Graduation	Tolerance	Shelf pack	Case
MICRO-BURET1S	1ml	0.005ml	±0.10ml		1
MICRO-BURET2S	2ml	0.01ml	±0.10ml		1
MICRO-BURET5S	5ml	0.02ml	±0.10ml		1
MICRO-BURET10S	10ml	0.02ml	±0.20ml		1

●Metal base stand for buret is not included. It is available separately.

Buret, Micro, Amber, without Stand

Code No	Specifications			Packing Quantity	
	Capacity	Graduation	Tolerance	Shelf pack	Case
B-MICRO-BURET1S	1ml	0.005ml	±0.10ml		1
B-MICRO-BURET2S	2ml	0.01ml	±0.10ml		1
B-MICRO-BURET5S	5ml	0.02ml	±0.10ml		1
B-MICRO-BURET10S	10ml	0.02ml	±0.20ml		1

●Metal base stand for buret is not included. It is available separately.



Metal base stand for Micro Buret

Code No	Description	Packing Quantity	
		Shelf pack	Case
BURET-DA1	Metal Base for Micro Buret		1