



Pipette, Measuring, Mohr Type

Code No.	Specifications				Packing Quantity	
	Capacity	Graduation	Tolerance	Length	Shelf pack	Case
M-0.1-PIPET	0.1ml	0.001ml	±0.005ml	350mm		10
M-0.2-PIPET	0.2ml	0.002ml	±0.005ml	350mm		10
M-0.5-PIPET	0.5ml	0.005ml	±0.01ml	350mm		10
M-1-PIPET	1ml	0.01ml	±0.015ml	360mm		10
M-2-PIPET	2ml	0.02ml	±0.02ml	360mm		10
M-5-PIPET	5ml	0.05ml	±0.03ml	375mm		10
M-10-PIPET	10ml	0.1ml	±0.05ml	385mm		10
M-25-PIPET	25ml	0.1ml	±0.1ml	455mm		1

●For pipetter spoits, see 7607 spoit bulb(P.33)



Pipette, Serological

Code No.	Specifications				Packing Quantity	
	Capacity	Graduation	Tolerance	Length	Shelf pack	Case
S-0.1-PIPET	0.1ml	0.001ml	±0.005ml	350mm		10
S-0.2-PIPET	0.2ml	0.002ml	±0.005ml	350mm		10
S-0.5-PIPET	0.5ml	0.005ml	±0.01ml	350mm		10
S-1-PIPET	1ml	0.01ml	±0.015ml	360mm		10
S-2-PIPET	2ml	0.02ml	±0.02ml	360mm		10
S-5-PIPET	5ml	0.05ml	±0.03ml	375mm		10
S-10-PIPET	10ml	0.1ml	±0.05ml	385mm		10
S-25-PIPET	25ml	0.1ml	±0.1ml	455mm		1

●For pipetter spoits, see 7607 spoit bulb(P.33)




Pipette, Volumetric, Transfer

Code No	Specifications			Packing Quantity	
	Capacity	Tolerance	Length	Shelf pack	Case
7100-0.5-EX	0.5ml	±0.005ml	300mm		10
7100-1-EX	1ml	±0.008ml	350mm		10
7100-2-EX	2ml	±0.01ml	350mm		10
7100-5-EX	5ml	±0.015ml	380mm		10
7100-10-EX	10ml	±0.02ml	440mm		10
7100-20-EX	20ml	±0.03ml	470mm		10
7100-25-EX	25ml	±0.03ml	480mm		10
7100-50-EX	50ml	±0.05ml	540mm		10



Pipette for Culture

Code No.	Specifications				Packing Quantity	
	Capacity	Graduation	Tolerance	Length	Shelf pack	Case
SC-PIPET1	1ml	0.01ml	±0.03ml	350mm		10
SC-PIPET2	2ml	0.02ml	±0.04ml	350mm		10
SC-PIPET5	5ml	0.05ml	±0.06ml	342mm		10
SC-PIPET10	10ml	0.1ml	±0.1ml	342mm		10
SC-PIPET20	20ml	0.2ml	±0.1ml	377mm		1






Volumetric Ware Pipette



7098 Pipette, Microliter, Lang-Levy

Code No	Specifications		Packing Quantity	
	Capacity	Tolerance	Shelf pack	Case
7098A-PIPET5	5 μ l	$\pm 1\%$	1	5
7098A-PIPET10	10 μ l	$\pm 1\%$	1	5
7098A-PIPET20	20 μ l	$\pm 1\%$	1	5
7098A-PIPET30	30 μ l	$\pm 1\%$	1	5
7098A-PIPET50	50 μ l	$\pm 1\%$	1	5
7098A-PIPET100	100 μ l	$\pm 1\%$	1	5
7098A-PIPET200	200 μ l	$\pm 1\%$	1	5
7098A-PIPET300	300 μ l	$\pm 1\%$	1	5
7098A-PIPET500	500 μ l	$\pm 1\%$	1	5



7607 Spoit Bulb, Silicon Rubber

Code No	Specifications		Packing Quantity	
	Bulb Capacity	Hole diameter	Shelf pack	Case
7607SPOIT1	1.5ml	5.5mm		10
7607SPOIT2	3.0ml	6.5mm		10
7607SPOIT3	4.5ml	6.5mm		10
7607SPOIT5	7.5ml	7.0mm		10
7607SPOIT10	15.0ml	8.0mm		10

- This spoit can be sterilized in an autoclave.
- Suitable for disposable pipette, SC pipette, etc.
- Suitable for 1~5ml M · S · 7100 pipettes.



7608 Spoit Bulb, Silicon Rubber

Code No	Specifications		Packing Quantity	
	Bulb Capacity	Hole diameter	Shelf pack	Case
7608SPOIT1	1.5ml	4.0mm		10
7608SPOIT2	3.0ml	4.0mm		10
7608SPOIT3	4.5ml	4.0mm		10
7608SPOIT5	7.5ml	4.0mm		10
7608SPOIT10	15.0ml	6.0mm		10

- This spoit can be sterilized in an autoclave.
- Suitable for M · S · 7100 pipettes.